

Insulomate



User Instructions

When u first buy the product.

Fill the syringe and the cannula part with insulin.

1

2

Open the device and place the syringe in the position shown.

Connect the cannula to the needle of the syringe.

3

4

Connect the butterfly cannula to the patient. (Obtain medical professionals assistance if necessary)

Make sure the end of the syringe is in contact with the slider. (Shown in diagram)

5

6

Close the device

Power on the device.

7

8

Confirm the device is switched on by the power light.

Wear the device on the body as required.

9

When u first buy the product.

Open the device and take the syringe out.

1

2

Refill the syringe.

Place the syringe so that the end is in contact with the slider as shown in image.

3

4

Close the device.

Press the reset button.

5

The device starts timing once you press the reset button. (Ex: First dose will be injected after 8 hours)g

In an Emergency

The background of the entire slide is a photograph of a medical device, likely a syringe pump. A clear plastic syringe is connected to a blue tube. The device has a control panel with several buttons and a small display screen. The image is slightly blurred and has a green tint.

Press the emergency button.

12 units will be injected once you press the button. The 8 hours start from the beginning once you press this button.

This can be pressed if you want a dose as soon as you refill the syringe.